

### Individual Package



SL901A



### Bulk Package



SL901A



### Product Application

Our disposable virus specimen collection tube are used for the collection and transportation of clinical, Influenza, Avian Influenza, Hand, Foot and Mouth Disease, Measles, Chlamydia, Mycoplasma, and Urea plasma specimens.

### Sample Requirements

1. Swab and Tube are sterilized using High Energy Electron Beam Irradiation, which can ensure the biological safety of the sampled person
2. Lot number located on each kit to allow for full specimen traceability Kits can be stored up to 2 year at room temperature
3. Expiration date on each tube, swab and kit Biohazard Bags with absorbent pad are available
4. Samples available to test and verify
5. Bulk packaging available to save you money

### Individual Package

50 kits/box

16 box/Carton

Box Size:

225\*145\*165mm

weight:0.7kg

Carton Size:

600\*480\*350mm

weight:12.5kg

### Bulk Package

50 kits/box

20 box/Carton

Box Size:


225\*120\*165mm

weight:0.65kg

Carton Size:

620\*470\*350mm

weight:14kg



RJ-1

Nasal Swab



Oral Swab



### Features

The sampling swab adopts polyamide (nylon) flocked, which can collect virus samples better than the traditional cotton swab head, and can elute samples more effectively, reduce the sample loss, and improve the detection sensitivity.

Product Name	Disposable Sampling Swab
Standard	Nasal Swab or Oral Swab
Swab Breakpoint	80mm,30mm
Swab Size	150 mm
Swab Materials	Swab tip:Flocked Nylon,Swab Stick:Food Grade ABS
Sterile Method	High Energy Electron Beam Irradiation
Certification	CE,NMPA,SGS,MSDS etc
OEM&ODM	Support
Package	450pcs/box,16boxes/Carton Carton Size:60*48*35cm